



SPECIAL PACKAGE



BENCHTOP METER

NEW

GET IT NOW

Benchtop multiparameter meter

Parameters that can be measured: Dissolved Oxygen, pH, mV, Temp, Conductivity, TDS, Salt, ORP (optional probe required)

Featured with magnetic stirrer and electrode suitable for OUR, SOUR and BOD testing

Featured with "Max and Min" mode, facilitating simplified testing procedures for OUR and SOUR tests.

Output data logging to PC directly via USB allowing users to visualize performance data in graphical form.

Memory data: 150 data

Switchable temperature selection of °C or °F

Easy replacement of electrode

REMARKS: OPTIONAL FOR PHOTOMETER IN BODSET AND BODCODSET

<https://csqanalytics.com>

info@csqanalytics.com

06-651 9118

CSQ[®]

Analytics

| Commitment | Satisfaction | Quality |