



# A combination of BOD & COD Set

**Item Code: BODCODset**

**This package incl. photometer PF-12PLUS, Heating Block VARIO Mini, Incubator, BOD & COD test kits, 10mL micropipette & BOD accessories**

## Programmable Heating & Refrigerated Incubator

<b>Volume</b>	50L
<b>Temperature range</b>	5-65°C
<b>Temperature resolution</b>	0.1°C
<b>Temperature fluctuation</b>	± 1°C
<b>Inner chamber size (D*W*H)</b>	50 x 350 x 400 mm

## Heating Block NANOCOLOR VARIO Mini

<b>Type</b>	Heating block for chemical and analytical digestion
<b>Number of bore holes</b>	6 (16 mm outer diameter)
<b>Display</b>	Graphic display 128 x 64 Pixel
<b>Temperatures</b>	6 preprogrammed temperatures (70°C, 100°C, 120°C, 148°C, 150°C, 160°C)
<b>Temperature Stability</b>	+/- 1°C (according to DIN, EN, ISO and EPA methods)

## Compact photometer PF-12Plus

<b>Type</b>	Filter photometer with microprocessor control, self-test & auto-calibration
<b>Light source</b>	Xenon Lamp
<b>Wavelengths</b>	345nm / 436nm / 470nm / 540nm / 585nm / 620nm / 690nm / 860nm
<b>Measuring Modes</b>	More than 100 preprogrammed test
<b>Cuvette Holders</b>	Round tubes 16 mm OD
<b>Memory</b>	1000 results, GLP conform
<b>Display</b>	Backlit graphic display, 128 x 64 pixel
<b>Operation</b>	Self-explanatory menu guidance, plastic foil keypad test selection via parameter lists
<b>Power Supply</b>	Via USB power supply, standard or rechargeable batteries or accupack
<b>Housing</b>	Shock-resistant, waterproof and dustproof, according to IP 68
<b>Dimensions</b>	215 X 100 X 65 mm
<b>Waranty</b>	2 years