

# BOD BASIC SET 2

REF: BODset2

Make your BOD5 testing easier than ever—  
all in one convenient package!



This package includes:

Benchtop Multimeter, Incubator, BOD5 nutrients  
& accessories and 10mL micropipette.



The images shown are for illustration purposes only & may not be an exact representation of the product



**BENCHTOP METER**



## BENCHTOP METER

1. Parameters that can be measured: pH, mV, Temp, Conductivity, TDS, Salt, Dissolved Oxygen, ORP (optional probe required)
2. Other Testing: OUR, SOUR, BOD testing
3. Featured with magnetic stirrer & electrode suitable for OUR and BOD testing
4. Featured with "Max & Min" mode, facilitating simplified testing procedures for OUR test.
5. Output data logging to PC directly via USB allowing user to visualize performance data in graphical form
6. Memory data: 150 data
7. Switchable temperature selection of °C or °F
8. Easy replacement of electrode



**PROGRAMMABLE HEATING & REFRIGERATED INCUBATOR**

## OSR PROGRAMMABLE HEATING & REFRIGERATED INCUBATOR

1. Volume : 50L
2. Voltage : AC 220V 50HZ or 110V 60HZ
3. Tem. range : 5-65°C
4. Tem. resolution : 0.1°C
5. Tem. fluctuation : ± 1°C
6. Power : 500W
7. Inner chamber size (D\*W\*H)mm : 350 x 350 x 400
8. Exterior size (D\*W\*H)mm : 560 x 460 x 910
9. Package size (D\*W\*H)mm : 610 x 530 x 1060
10. NW / GW (kg) : 45/53

✉ info@csqanalytics.com

🌐 <https://csqanalytics.com>

☎ 06-651 9118 / 017-928 3580 (WhatsApp)

**GET IT NOW!**

**CSQ**<sup>®</sup> Analytics  
| Commitment | Satisfaction | Quality |