

ADVANCED FILTRATION SOLUTIONS

MN615A

Types: **Qualitative filter paper**

Basis Weight (g/m²): **85**

Thickness (mm): **0.2**

Filtration Speed (s): **20**

Average Retention Capacity (µm): **4-12**

Binder: **None**

Equivalent to filter paper: **No.1**

Properties: **Medium fast filtration**

Applications: **General Filtration**

MN439

Types: **Glass fiber filter**

Basis Weight (g/m²): **66**

Thickness (mm): **0.42**

Filtration Speed (s): **13**

Max. Temperature (°C): **550**

Average Retention Capacity (µm): **1.5**

Binder: **None**

Equivalent to filter paper: **934-AH**

Properties: **Fast filtration**

Applications: **Total Suspended Solid Analysis**

***Experience
the Power of
Superior Filtration***

***Effortless
Performance,
Maximum Efficiency***

***Get Yours Today & Transform
Your Lab's Efficiency!***

Website: <https://csqanalytics.com>

Email: info@csqanalytics.com

Tel: 06-651 9118

CSQ[®] Analytics
| Commitment | Satisfaction | Quality |