


Analyze Smarter

ISOLAB Benchtop Meter

PERFECT for BOD testing!

- Other Testing: **DO, OUR, SOUR, pH, mV, Temp, Conductivity, TDS, Salt, Dissolved Oxygen, ORP** (optional probe required)
- Featured with magnetic stirrer and electrode suitable for BOD testing
- Output data logging to PC directly via USB allowing users to visualize performance data in graphical form.
- Memory data: 150 data
- Switchable temperature selection of °C or °F
- Easy replacement of electrode



Other items required for BOD testing



The images shown are for illustration purposes only & may not be an exact representation of the product

CSQ[®] | Analytics
| Commitment | Satisfaction | Quality |

info@csqanalytics.com

<https://csqanalytics.com>

06-651 9118 / 017-928 3580 (WhatsApp)